

Sample Treatment Parameters

Erbium-glass laser — reference material

Skin concern	Pulse width (ms)	Energy (W)	Spacing (mm)	Operating mode	Expected clinical endpoint
Pigmentation	1.0	8–12	0.7–1.0	QCW	Spot whitening of pigmented areas
Acne scars	1.0	10–14	0.7–1.0	QCW	Whitening of the skin, moderate edema
Enlarged pores	1.0	8–12	0.7	QCW	Mild edema and erythema
Fine wrinkles	1.0	8–12	0.7	QCW	Mild edema and erythema
Hair regeneration	0.5	4–6	0.7	QCW	Mild edema and erythema
Skin rejuvenation — mild	1.0	4–6	0.7–0.8	CW	Mild edema and erythema
Skin rejuvenation — moderate	1.5	6–8	0.7–0.8	CW	Mild edema and erythema
Skin rejuvenation — intensive	2.0	8–10	0.7–0.8	CW	Moderate edema and erythema
Active ingredient delivery	1.0	4–6	0.7–1.0	CW	Mild erythema

Note: The values shown are for reference only. Parameters must be selected according to the user manual for the specific model, the treatment type, the treatment area, and the patient's skin condition. The appropriate settings must be confirmed by a qualified professional before treatment.